



Material description: THYMOLPHTHALEIN

CAS No: 125-20-2

Supplied as: Product code Description  
CHE3782 THYMOLPHTHALEIN 5g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change Clear to Blue					pH 9.3-10.5

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1997/MS/19890101