

Material description: TOLUENE (glass distilled)

CAS No: 108-88-3

Supplied as: Product code Description
CHE3798 TOLUENE (glass distilled) 2.5L.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99.8 | | | |
| Water | % w/w | | 0.03 | | |
| Transmission; 284nm | % | 13 | | | |
| Transmission; 300nm | % | 78 | | | |
| Transmission; 310nm | % | 90 | | | |
| Transmission; 400nm | % | 99 | | | |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2937/MS/19890101