



Material description: TOLUIDINE BLUE O (C.I. 52040)

CAS No: 92-31-9

Supplied as: Product code Description  
CHE3806 TOLUIDINE BLUE O (C.I. 52040) 5g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wavelegth of Max Absorption	nm	628	633		
Loss on drying @110°C	%		8.00		

Revision date: 09/09/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3104/MS/20210909