

Material description: TRICHLOROACETIC ACID

CAS No: 76-03-9

Supplied as: Product code Description
CHE3810 TRICHLOROACETIC ACID 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.0			
Sulphate; as SO ₄	%		0.10		
Nitrate; as NO ₃	%		0.001		
Iron	%		0.0005		
Phosphate; as PO ₄	%		0.001		

Revision date: 08/02/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0710/MS/20210208