

Material description: TWEEN 80
CAS No: 9005-65-6

<i>Supplied as:</i>	Product code	Description
	CHE3854	TWEEN 80 1L.
	CHE5008	TWEEN 80 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wt. per ml.	g			1.1	
Water	% w/w		3		
Hydroxyl value	mg KOH/g	65	80		
Saponification value	mg KOH/g	45	55		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2790/MS/19890101