

Material description: UREA pure  
CAS No: 57-13-6

<i>Supplied as:</i>	Product code	Description
	CHE3866	UREA pure 500g.
	CHE3868	UREA pure 2Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98.5			
Biuret	%		1.0		
Moisture	%		0.5		

Revision date: 21/02/2012

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0751/MS/20120221