

Material description: WATER DEIONIZED

CAS No: 7732-18-5

<i>Supplied as:</i>	Product code	Description
	CHE3876	WATER DEIONIZED 25L.
	CHE3878	WATER DEIONIZED 5L.
	CHE3880	WATER DEIONIZED 2.5L.
	CHE4092	WATER DEIONIZED 200L.
	CHE4110	WATER DEIONIZED 10L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Conductivity @ 25°C	µS/cm		5.00		
Total Organic Carbon	ppb		500		

Revision date: 21/01/2016

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0757/MS/20160121