

Material description: WATER HPLC

CAS No: 7732-18-5

<i>Supplied as:</i>	Product code	Description
	CHE3882	WATER HPLC 2.5L.
	CHE4156	WATER HPLC 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Conductivity @ 25°C	µS/cm		1.00		
Total Organic Carbon	ppb		500		
Free acid	%		0.001		
Free alkali	%		0.001		
Non-volatiles	%		0.0002		
Abs. of peak @245 nm	AU		0.002		

Revision date: 07/02/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2674/MS/20200207