

Material description: WHITE SPIRIT
CAS No: 8052-41-3

<i>Supplied as:</i>	Product code	Description
	CHE3884	WHITE SPIRIT 1L.
	CHE3886	WHITE SPIRIT 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Appearance					Clear colourless, free from suspended material
Wt. per ml.	g			0.77	
Colour	Hazen		25		
Distillation range	°C	150	205		
Aromatics	%		25		

Revision date: 25/06/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1519/MS/20090625