

Material description: XYLENE A.R. (mixture of isomers)

CAS No: 1330-20-7

<i>Supplied as:</i>	Product code	Description
	CHE3898	XYLENE A.R. (mixture of isomers) 2.5L.
	CHE3900	XYLENE A.R. (mixture of isomers) 25L.
	CHE4048	XYLENE A.R. (mixture of isomers) 1L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	97.0			
Colour	Hazen		20.0		
Boiling range	°C	138	142		
Wt. per ml.	g	0.863	0.872		
Water	% w/w		0.05		

Revision date: 22/09/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0760/MS/20210922