

Material description: ZINC CHLORIDE ANHYDROUS pure

CAS No: 7646-85-7

<i>Supplied as:</i>	Product code	Description
	CHE3910	ZINC CHLORIDE ANHYDROUS pure 500g.
	CHE3912	ZINC CHLORIDE ANHYDROUS pure 3Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97.5			
Zinc oxide	%		2.20		
Copper	%		0.002		
Lead	%		0.005		

Revision date: 09/12/2015

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0762/MS/20151209