

Material description: ZINC NITRATE 6H<sub>2</sub>O pure  
CAS No: 10196-18-6

<i>Supplied as:</i>	Product code	Description
	CHE3924	ZINC NITRATE 6H <sub>2</sub> O pure 250g.
	CHE3926	ZINC NITRATE 6H <sub>2</sub> O pure 500g.

### Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Chloride	%		0.01		
Lead	%		0.005		

Revision date: 30/09/2015

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1523/MS/20150930