



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: ZINC SULPHATE 6H₂O pure

CAS No: 7446-19-7

<i>Supplied as:</i>	Product code	Description
	CHE3932	ZINC SULPHATE 6H ₂ O pure 500g.
	CHE3934	ZINC SULPHATE 6H ₂ O pure 1Kg.
	CHE3936	ZINC SULPHATE 6H ₂ O pure 25Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97.3			
Chloride	%		0.25		
Iron	%		0.01		

Revision date: 03/10/1997

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0766/MS/19971003