

Material description: TETRAHYDROFURAN HPLC  
CAS No: 109-99-9

*Supplied as:* Product code Description  
CHE3960 TETRAHYDROFURAN HPLC 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			
Water content	%		0.02		
UV Abs.: 400nm			0.0050		
UV Abs.: 350nm			0.0050		
UV Abs.: 300nm			0.0200		
UV Abs.: 250nm			0.1800		
UV Abs.: 212nm			1.0000		

Revision date: 22/03/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0675/MS/20170322