

Material description: POTASSIUM SULPHATE pure
CAS No: 7778-80-5

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|-------------------------------|
| | CHE3092 | POTASSIUM SULPHATE pure 500g. |
| | CHE3982 | POTASSIUM SULPHATE pure 100g. |
| | CHE3984 | POTASSIUM SULPHATE pure 1Kg. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 98 | | | |
| Loss on drying @ 130°C | % | | 1.00 | | |

Revision date: 19/01/2023

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0550/MS/20230119