

Material description: POTASSIUM SULPHATE pure  
CAS No: 7778-80-5

<i>Supplied as:</i>	Product code	Description
	CHE3092	POTASSIUM SULPHATE pure 500g.
	CHE3982	POTASSIUM SULPHATE pure 100g.
	CHE3984	POTASSIUM SULPHATE pure 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Loss on drying @ 130°C	%		1.00		

Revision date: 19/01/2023

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0550/MS/20230119