

Material description: POTASSIUM HYDROGEN PHTHALATE A.R.

CAS No: 877-24-7

<i>Supplied as:</i>	Product code	Description
	CHE2970	POTASSIUM HYDROGEN PHTHALATE A.R. 250g.
	CHE3994	POTASSIUM HYDROGEN PHTHALATE A.R. 100g.
	CHE3996	POTASSIUM HYDROGEN PHTHALATE A.R. 500g.
	CHE3998	POTASSIUM HYDROGEN PHTHALATE A.R. 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.95			
pH; 0.05M @ 25°C	pH	4.00	4.02		
Chloride	%		0.003		
Iron	%		0.0005		
Lead	%		0.0005		
Sodium	%		0.005		
Insoluble matter	%		0.005		
Sulphur compounds	%		0.002		

Revision date: 26/10/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2047/MS/20091026