

Material description: POTASSIUM HYDROGEN PHTHALATE A.R.

CAS No: 877-24-7

<i>Supplied as:</i>	Product code	Description
	CHE2970	POTASSIUM HYDROGEN PHTHALATE A.R. 250g.
	CHE3994	POTASSIUM HYDROGEN PHTHALATE A.R. 100g.
	CHE3996	POTASSIUM HYDROGEN PHTHALATE A.R. 500g.
	CHE3998	POTASSIUM HYDROGEN PHTHALATE A.R. 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.95			
pH; 0.05M @ 25°C	pH	4.00	4.02		
Chloride	%		0.003		
Iron	%		0.0005		
Lead	%		0.0005		
Sodium	%		0.005		
Insoluble matter	%		0.005		
Sulphur compounds	%		0.002		

Revision date: 26/10/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2047/MS/20091026