

Material description: PHENOLPHTHALEIN pure

CAS No: 77-09-8

<i>Supplied as:</i>	Product code	Description
	CHE2838	PHENOLPHTHALEIN pure 500g.
	CHE2840	PHENOLPHTHALEIN pure 250g.
	CHE2842	PHENOLPHTHALEIN pure 5g.
	CHE2844	PHENOLPHTHALEIN pure 100g.
	CHE4042	PHENOLPHTHALEIN pure 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Melting range	°C	258	263		
Colour change	pH				pH 8.3-10.0

Revision date: 09/05/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0508/MS/20180509