

Material description: XYLENE A.R. (mixture of isomers)

CAS No: 1330-20-7

<i>Supplied as:</i>	Product code	Description
	CHE3898	XYLENE A.R. (mixture of isomers) 2.5L.
	CHE3900	XYLENE A.R. (mixture of isomers) 25L.
	CHE4048	XYLENE A.R. (mixture of isomers) 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	97.0			
Colour	Hazen		20.0		
Boiling range	°C	138	142		
Wt. per ml.	g	0.863	0.872		
Water	% w/w		0.05		

Revision date: 22/09/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0760/MS/20210922