



Material description: PURECHLOR 1
CAS No: Mixture

Supplied as: Product code Description
CHE4052 PURECHLOR 1 1L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Density @20°C	g/ml	1.006	1.024		
Available chlorine	% w/v	1.0	1.1	1.0	

The above product is a synthetic mineral chemical. The raw material used in the manufacture is manufactured by our raw material supplier by absorbing Chlorine gas in in caustic soda. It is not derived from, nor does it contain, any products of animal origin. The raw material manufacturing plants are dedicated to the manufacture of Sodium hypochlorite and no animal materials are used when cleaning operations are conducted. The 1% Sodium hypochlorite solution is manufactured by us using the raw material detailed above and involves a simple dilution in water, within a manufacturing vessel. No animal products are used in any of our manufacturing vessels or our cleaning procedures.

Revision date: 04/02/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4809/MS/20200204