

Material description: SULPHURIC ACID 96% w/w pure
CAS No: 7664-93-9

Supplied as: Product code Description
CHE3942 SULPHURIC ACID 96% w/w pure 2.5L.
CHE4056 SULPHURIC ACID 96% w/w pure 10L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	95	97		
Residue on ignition	% w/w		0.01		
Free sulphur dioxide	ppm		20		
Iron	ppm		25		
Lead	ppm		1		
Mercury	ppm		0.1		

Revision date: 09/06/2022 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4125/MS/20220609