

Material description: OLEIC ACID pure
CAS No: 112-80-1

Supplied as: Product code Description
CHE5002 OLEIC ACID pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Free fatty acids ;as Oleic	%	97.5	102.5		
Acid value	mg KOH/g	194	204		
Iodine value	g I ₂ /100g	80	100		
Congeaing point	°C			10	

Revision date: 20/03/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0494/MS/20170320