



Material description: FORMALDEHYDE SOLUTION 40% w/w
CAS No: 50-00-0

Supplied as: Product code Description
CHE5032D FORMALDEHYDE SOLUTION 40% w/w 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	39	41		
Wt. per ml.	g			1.08	
Iron	%		0.002		
Acidity; formic	%		0.05		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2995/MS/19890101