



Material description: ZINC ACETATE 2H2O pure
CAS No: 5970-45-6

Supplied as: Product code Description
CHE5038 ZINC ACETATE 2H2O pure 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98	101		

Revision date: 01/03/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0761/MS/20170301