

Material description: 3M POTASSIUM CHLORIDE FILLING SOLUTION

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE2736	3M POTASSIUM CHLORIDE FILLING SOLUTION 100ml.
	CHE5062	3M POTASSIUM CHLORIDE FILLING SOLUTION 5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L			3.00	
Density @20°C	g/ml	1.130	1.134	1.132	

Revision date: 28/10/2019

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3124/MS/20191028