

Material description: XYLENE (packed in safety coated glass)

CAS No: 1330-20-7

*Supplied as:* Product code Description  
CHE5120 XYLENE (packed in safety coated glass) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range	°C	138	142		
Wt. per ml.	g	0.862	0.872		
Colour	Hazen		20		
Sulphur	%		0.0001		

Revision date: 03/10/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4615/MS/20171003