

Material description: XYLENE A.R. (packed in safety coated glass)

CAS No: 1330-20-7

Supplied as: Product code Description
CHE5122 XYLENE A.R. (packed in safety coated glass) 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range	°C	138	142		
Wt. per ml.	g	0.863	0.872		
Copper	%		0.00005		
Iron	%		0.0001		
Lead	%		0.00005		
Water	% w/w		0.05		
Non-volatiles	%		0.002		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4622/MS/19890101