



Material description: NITRIC ACID 1% v/v

CAS No: 7697-37-2

Supplied as:	Product code	Description
	CHE6018	NITRIC ACID 1% v/v 1L.
	CHE6018	NITRIC ACID 1% v/v 5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.232	0.238		

Revision date: 14/11/2019

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0373/MS/20191114