

Material description: SODIUM SILICATE POWDER pure
CAS No: 1344-09-8

| | | |
|---------------------|--------------|-----------------------------------|
| <i>Supplied as:</i> | Product code | Description |
| | CHE3530 | SODIUM SILICATE POWDER pure 500g. |
| | CHE6019 | SODIUM SILICATE POWDER pure 25Kg. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-------------------|-------|---------------|---------------|---------------|-------|
| SiO ₂ | % | | | 53 | |
| Na ₂ O | % | | | 26 | |

Revision date: 21/02/1991 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1428/MS/19910221