



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: BARIUM SULPHATE pure

CAS No: 7727-43-7

| Supplied as: | Product code | Description                |
|--------------|--------------|----------------------------|
|              | CHE6044      | BARIUM SULPHATE pure 3Kg.  |
|              | CHE6046      | BARIUM SULPHATE pure 500g. |

## Specification

| Parameter                  | Units | Minimum value | Maximum value | Typical value | Notes   |
|----------------------------|-------|---------------|---------------|---------------|---------|
| Acid insoluble matter      | %     |               | 0.3           |               |         |
| Soluble Ba compounds       | %     |               | 0.001         |               |         |
| Acidity/alkalinity; Litmus |       |               |               |               | Neutral |

Revision date: 18/09/1992

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 4105/MS/19920918