



Key Features

- Most Energy Efficient
- Low Noise
- Reliable Performance

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l578j.html

Specifications

Technical Data

Type	-86°C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	HC
Noise Level (dB)	50

Performance

Cooling performance(°C)	-86
Temperature Range(°C)	-40~-86

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
---------------------------	------------

Power(W)	900
Electrical Current(A)	9
Dimensions	
Capacity(L/Cu.Ft)	578/20.4
Net/Gross Weight(approx)	300/330(kg) 661.4/727.5(lbs)
Interior Dimensions(W*D*H)	620*716*1310(mm) 24.4*28.2*51.6(in)
Exterior Dimensions(W*D*H)	895*980*1960(mm) 35.4*38.6*78.0(in)
Packing Dimensions(W*D*H)	950*1055*2125(mm) 37.4*41.5*83.7(in)
Container load (20'/40'/40'H)	12/24/24
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	N/A
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	Y
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
-------------------	-------------------------	------------------------------------	----------------------	----------------------	--

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate

CE
