



Key Features

- Special designed Slot For Sign Card
- Adjustable Shelves
- Five Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/pharmacy_refrigerator/hyc-1378.html

Specifications

Technical Data

Type Pharmacy Refrigerator

Cabinet Type Upright

Climate Class N

Cooling Type Forced air cooling

Defrost Mode Auto

Refrigerant CFC-Free

Noise Level (dB) 43

Performance

Temperature Range(°C) 2~8

Control

Controller Microprocessor

Display LED

Electrical Data

Power Supply(V/Hz) 220-240/50/208-230/60/115/60

Power(W) 800 800 900

Electrical Current(A)	4.2 4.0 10
------------------------------	------------

Dimensions

Capacity(L/Cu.Ft)	1378/48.7
--------------------------	-----------

Net/Gross Weight(approx)	310/360(kg) 683.4/793.7(lbs)
---------------------------------	---------------------------------

Interior Dimensions(W*D*H)	1320*700*1500(mm) 52.0*27.6*59.1(in)
-----------------------------------	---

Exterior Dimensions(W*D*H)	1440*875*1980(mm) 56.7*34.4*78.0(in)
-----------------------------------	---

Packing Dimensions(W*D*H)	1600*1010*2100(mm) 63.0*39.8*82.7(in)
----------------------------------	--

Container load (20'/40'/40'H)	6/12/12
--------------------------------------	---------

Functions

High/Low Temperature	Y
-----------------------------	---

Remote Alarm	Y
---------------------	---

Power Failure	Y
----------------------	---

Sensor Error	Y
---------------------	---

Low Battery	Y
--------------------	---

Door Ajar	Y
------------------	---

Accessories

Caster	Y
---------------	---

Foot	Y
-------------	---

Test Hole	Y
------------------	---

Shelves/Inner Doors	12/-
----------------------------	------

USB Interface	Y
----------------------	---

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate	CE
--------------------	----