



DATA SHEET

Product Code SLS1324 – SLS Select 2.2mL square well 'V' bottom plate with notches

About SLS1324

The 2.2-mL square well 'V' bottom plate has been developed to enable the base of the well to be in contact with heater-shaker blocks which improves the performance of the process. It has a well depth of 41.3mm and a total volume of 2.4mL.

It has the added benefit of notches on the sidewall for compatibility with KingFisher™ platforms.

The plate is manufactured in SLS Select's class 8 cleanrooms using a medical-grade polypropylene. This polypropylene helps achieve very low concentrations of extractable elements and has a temperature capability range of -80°C to +121°C.

This plate has been developed precisely to conform to the footprint dimensions defined by SLAS/ANSI standards. This guarantees its compatibility with both manual and automated sample handling systems, readers, and microplate washers.



United Kingdom:
slsinfo@scientific-labs.com

International:
slsexport@scientific-labs.com

Northern Ireland:
NI.sales@scientific-labs.com

Republic of Ireland:
ROI.sales@scientific-labs.com

 **SCIENTIFIC
LABORATORY
SUPPLIES**

