



# Product Data Sheet

## SLS5008

### SLS Select Filter Tips

20µl (1-20µl) Grad Barrier



Catalog #	Description
SLS5008	Natural, Polypropylene, Sterile, Filtered, Low Binding

Packaging Type	Packaging Details
Racked Tips	MLT GRD HNG FRAME 20 STDZN ST Ir 4800/C 960/P

#### Attributes



**Material:** Polypropylene  
**VolumeRange:** 10-20µl  
**Grad Lines:** 10

Certified RNase/DNase Free  
Certified Human DNA Free  
Certified PCR Inhibitor Free  
Certified Pyrogen Free