



Product Data Sheet

SLS8114

General Descriptions

Product Name	SLS Select centrifuge tube, 50mL
Further Description	Self-standing, skirted conical base Transparent Flat screw cap Sterile Manufactured from virgin ultra clear polypropylene to ensure good visibility Leak-resistant Maximum accuracy of graduations at $\pm 2\%$ White printed graduations with large white marking strip for labelling
Container Material	Polypropylene
Cap Material	HDPE
Packaging	Bagged

General Features

Graduations	White graduation lines up to 50mL
Sterility	ETO sterile (SAL = 10^{-6})
Tube Dims, OD x h	30 x 116mm
Satisfies the Requirement	IATA, -95kPa (14psi)
Product Category	IVD
DNase/RNase	Free
Certificates Certified	CE (98/79/EC) / FDA
Non-Pyrogenic	<0.125EU/mL
Autoclavable	Max. 121°C
Leak-Proof at	-95kPa (14psi)
Max RCF	9,500 Refers to liquids with a specific gravity of 1.0, using appropriate rotors, cushions, and adapters

