

# Cole-Parmer®

## Cole-Parmer® Stirrer Multiposition, ST-200-MP

- Three independently controlled stirring positions
- Powerful magnets for strong magnetic coupling
- Stainless steel top plate





## Cole-Parmer® Stirrer Multiposition, ST-200-MP

A space saving and economical unit with three independent stirring positions in a footprint only 600 x 270mm. Easily accommodates 3 x 2 litre beakers. The stainless steel top plate is robust and, unlike aluminium, does not produce eddy currents and so ensures a very powerful coupling and stirring action. Powerful magnets and motor give stirring speed up to 1500rpm and volumes up to 15 litres \*. Supplied complete with 3 x 25mm PTFE coated stir bars.

### Technical Specification

Specification	ST-200-MP
Plate material	Stainless steel
Plate dimensions, mm	550 x 210
Speed range, rpm	100 to 1500
Dimensions, mm (w x d x h)	606 x 268 x 94
Net weight, kg	7
Electrical supply	230V, 50Hz, 150W
IP Rating	31

### Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Stirrer, magnetic, 3 position and Allen key	51400-02	ST-200	ST-200-MP	SB161-3
Retort rod bracket	04805-96	ST-200	N/A	SR3
Retort rod, 600mm x 12mm diameter.	04805-91	ST-200	N/A	SR1

[coleparmer.com](http://coleparmer.com)

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials

an antylia scientific company

Find out more at [antylia.com](http://antylia.com)