

# Cole-Parmer®

## Cole-Parmer® Hotplate Stirrer, Multiposition, SHP-200

- Three independently controlled stirring/heating positions
- Separate “Hot” warning lights for each plate
- Powerful magnets for strong magnetic coupling





## Cole-Parmer® Hotplate Stirrer, Multiposition, SHP-200

A space saving and economical unit with 3 independent heating/stirring positions in a footprint only 600 x 270mm. Unit only requires one power point and is ideal for quality control applications where multiple samples require simultaneous heating and stirring, under the same conditions.

A “Hot” warning light for each plate will flash whenever its temperature is above 50°C and will operate even when the hotplate is turned off and connected to the mains. The top plates are a robust aluminium/ silicon alloy, providing even plate temperature. Powerful magnets and motor give stirring speed up to 1500rpm and volumes up to 15 litres.

\* Supplied complete with 3 x 25mm PTFE coated stir bars.

### Technical Specification

Specification	SHP-200
Plate material	Al/Si alloy
Individual plate dimensions, mm	160 x 160
Heater power, W	3 x 700
Maximum plate temperature, °C	325
Speed range, rpm	100 - 1500
Dimensions, mm (w x d x h)	600 x 270 x 110
Net weight, kg	11
Electricity supply	230V, 50Hz, 2250W
IP Rating	32

\* Based on water contained in a 20 litre glass bottom boiling flask.

### Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Stirrer, magnetic, 3 position with heating	51600-08	SHP-200	SHP-200-MP	SB162-3
Retort rod bracket	04805-96	SHP-200	N/A	SR3
Retort rod, 600mm x 12mm diameter.	04805-91	SHP-200	N/A	SR1

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials

an antylia scientific company

Find out more at [antylia.com](http://antylia.com)