

Cole-Parmer®

Cole-Parmer® Portable Magnetic Stirrer, ST-200

- Battery powered or mains supply
- Tough, compact and light weight
- Variable speed up to 1300rpm





Cole-Parmer® Portable Magnetic Stirrer, ST-200

Cole-Parmer® ST-200 portable, magnetic stirrer, powered by long life batteries or mains power this rugged little unit can be used almost anywhere. Ideal for use inside incubators, glove boxes or in the field.

Housed in a tough, chemically resistant ABS case and supplied complete with batteries. Operating indicator will dim when batteries run low.

Technical Specification

Specification	ST-200
Max speed, rpm	1300
Capacity, litres H ₂ O	1.5
Dimensions, mm (w x d x h)	150 x 160 x 70
Net weight, kg	1.3
Electricity supply	4 x alkaline batteries or 230V, 50Hz
Battery type	1.5V, size D, IEC No. LR20
Battery life	600 hours continuous use
IP rating	41

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Portable magnetic stirrer complete with 4 x LR20 batteries and 25mm PTFE stir bar	04801-05	ST-200	ST-200-P	SM27
Mains adapter, 9V a.c. for 230V, 50Hz supply with UK plug	04801-69	ST-200	N/A	SM27/1
Mains adapter, 9V a.c. for 230V, 50Hz supply with European 2 pin plug	04801-71	ST-200	N/A	SM27/2
Car battery adapter, 12V, Fitted with standard car cigarette/cigar plug	04801-73	ST-200	N/A	SM27/3

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials

an antylia scientific company

Find out more at antylia.com