




designed for scientists



R 1401 Propeller

/// Data Sheet

The R 1401 propeller head is designed for the use with the overhead stirrers EUROSTAR 20 high speed digital / control. It is ideally suited for applications which work with low to medium mixing viscosities that require axial flow pattern at high speed with minimal temperature increase.

Working range: 1 - 30 l

Rotor diameter: 55 mm

Note: R 1401 needs to be attached to the precision shaft R 6000 in order to operate. The shaft is not included in the



designed for scientists

scope of delivery of R 1401.



designed for scientists

Technical Data

Stirrer diameter [mm]	55
Weight [kg]	0.0357
Permissible ambient temperature [°C]	5 - 40