


IKA

designed for scientists



R 1827 Plate stand

/// Data Sheet

With slip resistant foil.

Diameter of support rod: 16 mm

Max. load: 5 kg

Dimensions (W x D): 200 x 316 mm

Height: 1000 mm

www.ika.com

Subject to technical changes





designed for scientists

Technical Data

Diameter [mm]	16
Dynamic load [kg]	5
Dimensions (W x H x D) [mm]	200 x 1000 x 316
Weight [kg]	6.326
Permissible ambient temperature [°C]	5 - 40