



### R 1300 Dissolver stirrer

Radial flow, for drawing the material to be mixed from the top and the bottom. High turbulence, high shearing forces. Particle reduction. Used at medium to high speeds.

Stirrer Ø: 80 mm

Shaft Ø: 8 mm

Shaft length: 350 mm

#### Technical Data

Stirrer diameter [mm]	80
Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
<b>Ident. No.</b>	<b>0000513500</b>