



R 301 Stirring shaft protection

/// Data Sheet

Prevents potential injuries at the rotating shafts and stirring elements. It can be directly attached to the stirring motors RW 20 digital, RW 28 digital and the EUROSTAR series overhead stirrer.

R 301.1 support holder is required to fix the R 301 stirring shaft protection when use with NANOSTAR 7.5 digital, MINISTAR, MICROSTAR series and Eurostar 400 overhead stirrer.

Length adjustment: 190 - 310 mm

Material: plexiglass



designed for scientists

Technical Data

Weight [kg]	0.204
Permissible ambient temperature [°C]	5 - 40