


IKA

designed for scientists



R 1376 Paddle stirrer

/// Data Sheet

Tangential flow, minimum turbulence, good heat exchange, gentle treatment of product. Used at low to medium speeds.

Stirrer Ø: 150 mm

Shaft Ø: 10 mm

Shaft length: 550 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	150
Shaft diameter [mm]	10
Shaft length [mm]	550
Speed max. [rpm]	800
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.6385
Permissible ambient temperature [°C]	5 - 40