


# 22s gauge, Small Hub RN Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 7758-03)



Small Hub Removable (RN) needles can be used with removable (RN)-style syringes that are 100  $\mu$ L and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

<b>Gauge</b>	22s gauge
<b>Point Style</b>	2
<b>Needle Length</b>	2 inch (51 mm)
<b>Autoclavable</b>	Yes
<b>Sales Quantity</b>	6 needles per pack
<b>Hub</b>	Small RN

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/needles/7758-03>