

# Lovibond® Water Testing Tintometer® Group



## Thermoreactor RD 125

For the Tube test digestion of



Part Number: 2418940

### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

### Technical Data

<b>Suitable Vials</b>	Round Cuvettes 16 mm
<b>Power Supply</b>	230 V / 50-60 Hz or 115 V / 50-60 Hz
<b>Auto – OFF</b>	Yes
<b>Beeper</b>	max. 88 dB (piezo buzzer)
<b>Suggestion Temperature</b>	100 / 120 / 150 °C
<b>Heating-up Rate</b>	20 °C - 150 °C
<b>Selected Time</b>	30 / 60 / 120 / 分钟和持续运行
<b>Environmental Conditions</b>	10 – 40 °C max. 85 % rel. humidity
<b>Compliance</b>	CE
<b>Dimensions</b>	219 x 171 x 248 mm
<b>Weight</b>	3.9 kg

### Accessory

<b>Title</b>	<b>Part Number</b>
Cuvette stand for 10 round cuvettes Ø 16 mm	418957
Pipette tips*, 1-5 ml (white) 100 pc.	419066
Pipette tips**, 0,1-1 ml (white), 1000 pc.	419073
Automatic pipette*, 1-5 ml	419076
Automatic pipette**, 0,1-1 ml	419077