

SLS Lab Pro 20T2 PurA-Q2

Complete All-in-One Type 2 Water System

- 1-15MΩ-cm Type 2 water
- Space saving footprint with integral 35L tank
- Make up rate 20L/hr
- Low energy consumption <40W
- Quiet internal boost pump and patented anti-vibration mountings
- Mercury-free option available
- First consumables included with unit
- Simple cartridge replacement



Code	Description
WAT2408	SLS Lab Pro 20T2 (PurA-Q2)
WAT2418	SLS Lab Pro 20T2 (PurA-Q2) with UV kit

Accessories

Code	Description
WAT2456	Wall mounting kit

Consumables

Code	Description
WAT2436	SLS Lab Pro pre-treatment module (RO) 20L/hr – *annual change
WAT2438	SLS Lab Pro DI module – *annual change
WAT2440	SLS Lab Pro air vent filter (tank) – *annual change
WAT2442	SLS Lab Pro UV lamp consumable for WAT2418 – *annual change

*Estimate only, depends on water usage and quality of feed water.

Model Specifications

Model	RO Make-up Rate	Inorganics Rejection Rate	Organics Rejection Rate
SLS Lab Pro 20T2	20L/Hr	1-15MΩ-cm	>99%
Bacteria Rejection Rate	pH	Daily Usage (Max)	
99%	Neutral	120L/day	



Applications

Type 2 Water Applications

Buffer & Media Preparations
 Sample Dilution & Reagent Preparation
 Spectrophotometry
 Protein Electrophoresis
 Cytology & Histology
 Glassware Washing & Rinsing

Dimensions			Feed Water Requirements		Services Required
Width	Height	Depth	Pressure	- 5 Bar	
485mm	500mm	330mm	Resistivity	<1400 μS-cm	
Weight (Empty)			Colloid Index	<3 SDI	
20kg			Free-chlorine & Iron	0.1 Mg/L	

00160